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Instruments and Software

· ChroZen GC/MS System

Item	Description	Part No.
Oven	ChroZen GC Mainframe Assembly for Mass Spectrometer	6701012500
Inlet	Capillary Inlet Assembly for ChroZen GC	6701012550
Detector	ChroZen MS for ChroZen GC incl. built-in turbomolecular pump - Single Quadropole - EI source (UEIS:Ultimated Efficiency Ion Source) - Turbo pump(240 L/sec) - Include fore pump and spares kit - Higher Sensitivity (S/N for OFN 2500:1)	6901012110
ChroZen PAL RSI Autosampler	PAL RSI 850 system without liquid or headspace option consisting of: 1 pc PALbase 850 X/Y unit with PAL Control Board 1 pc PAL RSI-zhead 1 pc PALterminal 1 pc Safety guard 1 pc PALpower and cables 1 pc Gas Purge Line for Headspace Technique built-in (Requires Mounting Kit for GC Model, additionally Tools, Syringes, Tray Holders or Stack and Racks)	PAL3-RSI-GH-S-9- YL2
	Mounting Kit for ChroZen GC	PAL3-Kit-YL6700
	Agitator for the incubation and agitation of up to 6x 20mL vials requires firmware version 1.3.0 or higher	PAL3-Agitator
ChroZen PAL RSI Optional Module (SPME Arrow)	Smart SPME Arrow Kit, consisting of: 1 pc SPME Arrow Tool 1 pc SMART SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix	PAL3-SPME-SArr-Kit
	SPME Arrow Conditioning Module for the conditioning of SPME arrows and SPME fibers Includes a second position for conditioning a replacement SPME arrow or SPME fiber requires firmware version 2.3 or higher	PAL3-SPME- ArrowCond
	Heatex Stirrer for intensive heating and stirring of 1x 20mL vial without stir bar Specific Heatex Stirrer insert for 10mL vial available (not included) requires firmware version 2.3 or higher (use in combination with existing PAL3-Agitator)	PAL3-HeatexStirrer
	Adaptation for split/splitless injector of YL6500 GC	Call
	SPME Arrow Liner ID 1.3mm, Pk of 3 pcs for SSL Injector of AGILENT GC 7890, only for SMPE Arrow with OD 1.1mm	ARRLIN13-GC7890-3
CDS	YL-Clarity software for single instrument of YL GC	5301011020
	MS module of YL-Clarity(Library 별도)	5301011180
	Autosampler control of YL-Clarity	5301011040
	Library(NIST/EPA/NIH 2017 edition)-306,622 spectra	NIST17-MS-LIB
Column	YL-5MS (30m, 0.25 mm, 0.25 μm)	1256120170

<p>Install. Option</p>	<p>Start-up kit includes (Without GC Capillary Column)</p> <ol style="list-style-type: none"> 1) Nuts and ferrules 1/8" 2) Nuts, 1/4" 3) Vespel ferrule, 1/4" 4) Union Tee, 1/8" 5) Septa, 11 mm, 50/pk 6) Tubing cutter 7) Monkey spanner 8) Wrench 1/2" & 9/16" 9) Wrench 1/2" & 7/16" 10) Wrench 3/8" & 7/16" 11) Wrench 1/4" & 5/16" 12) (+) screw driver 6x100 13) (+) screw driver 5x100 14) (-) screw driver 6x100 15) Leak detection fluid 16) 10 µL syringe 17) Copper tubing, 1/8" 	<p>1601011110</p>
<p>ACC</p>	<p>Big Universal Trap, 1/8" fttgs, Helium</p> <p>Sample Rack for 15 vials of 10/20mL, pack of 3 pcs Vials are not included</p> <p>Tray Holder for MTP und DW Plates, VT54, VT15 (CTC)</p> <p>YL PAL System Vial 20CV, 20ml Clear Glass with Label, designed for the YL PAL Autosampler. 75.5x22.5mm, 1st Class Hydrolytic Glass with flat finish for better sealing, fits ND18 Screw Caps, Pk of 100 Pcs</p> <p>YL PAL System Screw Cap 10CV and 20CV, designed for the YL PAL Autosampler. ND18 Magnetic, Silver, Silicone white/PTFE blue Septa 1.5mm, Pk of 100 pcs</p>	<p>RMSH-2</p> <p>PAL3-Kit-VT15</p> <p>PAL3-Trayplate</p> <p>Vial-20-ND18-CG-100</p> <p>Cap-ND18-St-SP15-100</p>

Reagents and Standards

Geosmin, 2-MIB Odor Compounds Mix
(200 µg/mL each component in
methanol)

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Table 1. GC/MS Condition

SPME System Conditions
Preconditioning time: 10min
Incubation time: 2min
Sample Extraction time: 30min
Extraction Temp.: 60°C
Sample Desorb time: 5min

Table 2. SPME System Condition

	Compound	Formula	Molecular Weight	Selected Ions, m/z
1	2-Methylisoborneol	C ₁₁ H ₂₀ O	168.28	95
2	Geosmin	C ₁₂ H ₂₂ O	182.31	112

Table 3. GC/MS Chromatographic Conditions


Fig 2. 2-MIB, Geosmin Concentration Chromatogram Overlay


Fig 3. 2-MIB, Geosmin 10ppt Chromatogram


Fig 4. 2-MIB Spectra


Fig 5. Geosmin Spectra



Fig 6. 2-MIB, Geosmin Calibration (1, 5, 10, 50,100ppt)

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