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Instruments and Software

ChroZen GC/MS System

Item	Description	Part No.
Oven	ChroZen GC Mainframe Assembly for Mass Spectrometer	6701012500
Inlet	Capillary Inlet Assembly for ChroZen GC	6701012550
Detector	ChroZen MS for ChroZenGC incl. built-in turbomolecular pump - Single Quadropole - EI source (UEIS:Ultimate Efficiency Ion Source) - Turbo pump(240 L/sec) - Include fore pump and spares kit - Higher Sensitivity (S/N for OFN 2500:1)	6901012110
Install. Option	Start-up kit (Without GC Capillary Column)	1601011110
YCM Liquid Autosampler	ChroZen PAL LSI system for liquid injection consisting of: 1 pc base 850 X/Y unit with PAL Control Board 1 pc head 1 pc Safety guard 1 pc power and cables 1 pc Liquid Syringe Tool 1 pc Standard Wash Module 2 pc Smart Liquid Syringe 10uL, Gauge 26S, conical, metal plunger 1 pc Trayplate 3 pc Sample Racks for 54 vials of 2 mL	6501011590
	Mounting Kit for ChroZen GC	PAL3-Kit-YI6700
CDS	YL-Clarity software for single instrument of YL GC	5301011020
	MS module of YL-Clarity(Library 별도)	5301011180
	Autosampler control of YL-Clarity	5301011040
	Library(NIST/EPA/NIH 2017 edition)-306,622 spectra	NIST17-MS-LIB
Column	YL 5ms (30m x 0.25mm i.d. x 0.25um)	1256120170
ACC	ChroZen PAL System Vial 2CV, 1.5mL Clear Glass with Label, designed for the YL PAL Autosampler. 12x32mm, 1st Class Hydrolytic Glass, fits ND9 Screw Caps, Pk of 100 Pcs	Vial-1.5-ND9-CG-100
	ChroZen PAL System Screw Cap 2CV, designed for the YL PAL Autosampler. ND9, Magnetic, Gold, Silicone/PTFE Septa 1.0mm, Pk of 100 pcs	Cap-ND9-St-SP10-100

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Table 2. GC/MS Conditions

	Compound	Molecular Weight	Selected Ions, m/z
1	Octylphenol	206	135, 107
2	4-n-Nonylphenol (SS)	220	107, 220
3	Phenanthrene-d10 (IS)	188	188, 187

Table 3. GC/MS Chromatographic Conditions


Fig 2. Chromatogram


Fig 3. Verification of Calibration Curve

R.T (min)	방법검출한계 (µg /L)	정량한계 (µg /L)	Accuracy (%)	Precision
16.061	0.1589	0.5062	98.68	0.05

Table 4. Validity of Test Method

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Reference

- EPA Method 8270E_ Semivolatile Organic Compounds By Gas Chromatography/Mass Spectrometry
- ASTM D7065-06_ Standard Test Method for Determination of Nonylphenol, Bisphenol A, p-tert-Octylphenol, Nonylphenol Monoethoxylate and Nonylphenol Diethoxylate in Environmental Waters by Gas Chromatography Mass Spectrometry
- 수질오염시험기준 ES 04614.1_옥틸페놀 기체크로마토그래피-질량분석법
- 환경시험 검사 QA/QC 핸드북_국립환경과학원(2011)

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